

Product datasheet

Anti-ZNHIT1 antibody ab68467

1 Image

Overview

Product name	Anti-ZNHIT1 antibody
Description	Mouse polyclonal to ZNHIT1
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	ZNHIT1 (NP_006340, 1 a.a. ~ 154 a.a) full-length human protein.
Positive control	Transfected 293T cell line.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 1X PBS, pH 7.2
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab68467** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		1/500 - 1/1000. Detects a band of approximately 18 kDa (predicted molecular weight: 18 kDa).

Target

Relevance

ZNHIT1, Zinc finger HIT domain containing protein 1, appears to play a role in p53-mediated apoptosis induction.

Images



All lanes : Anti-ZNHIT1 antibody (ab68467) at 1/500 dilution

Lane 1 : ZNHIT1 transfected 293 lysate

Lane 2 : Non-transfected 293 lysate

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP at 1/2500 dilution

Predicted band size: 18 kDa

Observed band size: 18 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors