

Product datasheet

Anti-TTC27 antibody ab68740

2 Images

Overview

Product name	Anti-TTC27 antibody
Description	Mouse polyclonal to TTC27
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Full length human protein (entry NP_060205, aa 1-843)
Positive control	Human pancreas; TTC27-transfected 293T cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 1X PBS, pH 7.2
Purity	Protein G purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab68740** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Predicted molecular weight: 97 kDa.

Target


Sequence similarities

Belongs to the TTC27 family.

Contains 3 HAT repeats.

Contains 7 TPR repeats.

Images




Anti-TTC27 antibody (ab68740) at 1/500 dilution + human pancreas tissue lysate at 25 μ g

Secondary

Goat Anti-Mouse IgG (H&L)-HRP conjugated at 1/2500 dilution

Predicted band size: 97 kDa

Observed band size: 88 kDa



All lanes : Anti-TTC27 antibody (ab68740) at 1/500 dilution

Lane 1 : TTC27 transfected lysate

Lane 2 : non transfected lysate

Lysates/proteins at 25 μ g per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP conjugated at 1/2500 dilution

Predicted band size: 97 kDa

Observed band size: 88 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors