

Product datasheet

Anti-TTC27 antibody ab68740

2 Images

Overview

<b>Product name</b>	Anti-TTC27 antibody
<b>Description</b>	Mouse polyclonal to TTC27
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	Full length human protein (entry NP_060205, aa 1-843)
<b>Positive control</b>	Human pancreas; TTC27-transfected 293T cell lysate

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Storage buffer</b>	Preservative: None Constituents: 1X PBS, pH 7.2
<b>Purity</b>	Protein G purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab68740** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		1/500 - 1/1000. Predicted molecular weight: 97 kDa.

Target

**Sequence similarities** Belongs to the TTC27 family.

Contains 3 HAT repeats.  
Contains 7 TPR repeats.

## Images




Anti-TTC27 antibody (ab68740) at 1/500 dilution + human pancreas tissue lysate at 25  $\mu$ g

### Secondary

Goat Anti-Mouse IgG (H&L)-HRP conjugated at 1/2500 dilution

**Predicted band size** : 97 kDa

**Observed band size** : 88 kDa



**All lanes** : Anti-TTC27 antibody (ab68740) at 1/500 dilution

**Lane 1** : TTC27 transfected lysate

**Lane 2** : non transfected lysate

Lysates/proteins at 25  $\mu$ g per lane.

### Secondary

Goat Anti-Mouse IgG (H&L)-HRP conjugated at 1/2500 dilution

**Predicted band size** : 97 kDa

**Observed band size** : 88 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <http://www.abcam.com/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors