

Product datasheet

Anti-DHDH antibody ab69435

[2 Images](#)

Overview

Product name	Anti-DHDH antibody
Description	Mouse polyclonal to DHDH
Host species	Mouse
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Full length human DHDH (NP_055290.1).
Positive control	Human colon tissue lysate and DHDH transfected 293T cell lysate.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 1X PBS, pH 7.2
Purity	Protein G purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab69435** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/1000. Detects a band of approximately 39 kDa (predicted molecular weight: 37 kDa).


Target

Tissue specificity

Small intestine.

Sequence similarities

Belongs to the Gfo/l dh/MocA family.

Images


Anti-DHDH antibody (ab69435) at 1/500 dilution + human colon tissue lysate at 25 µg

Secondary

Goat Anti-Mouse IgG (H&L)-HRP at 1/2500 dilution

Predicted band size: 37 kDa

Observed band size: 43 kDa



All lanes : Anti-DHDH antibody (ab69435) at 1/500 dilution

Lane 1 : DHDH transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP at 1/2500 dilution

Predicted band size: 37 kDa

Observed band size: 39 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors