

Product datasheet

Rat monoclonal SB87 α Anti-Rabbit IgG+IgM H&L preadsorbed ab80301

Overview

Product name	Rat monoclonal SB87 α Anti-Rabbit IgG+IgM H&L preadsorbed	
Target species	Rabbit	
Specificity	Reacts with Rabbit IgG and Rabbit IgM as demonstrated by ELISA; may also react with Ig from other species.	
Tested applications	Suitable for: Flow Cyt, ELISA	
Minimal cross-reactivity	Dog, Goat, Guinea pig, Hamster, Human, Monkey, Mouse, Pig, Rat, Sheep	more details
Immunogen	The details of the immunogen for this antibody are not available.	

Properties

Form	Liquid	
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.	
Storage buffer	pH: 8.20 Constituent: 100% Borate buffered saline	
	No preservatives or amine-containing buffer salts added.	
Purity	IgG fraction	
Clonality	Monoclonal	
Clone number	SB87 α	
Isotype	IgG2a	

Applications

Our [Abpromise guarantee](#) covers the use of **ab80301** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent dilution.

Application	Abreviews	Notes
ELISA		1/4000 - 1/8000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors