

## Product datasheet

# Rat monoclonal SB87a Anti-Rabbit IgG+IgM H&L preadsorbed ab80301

### Overview

<b>Product name</b>	Rat monoclonal SB87a Anti-Rabbit IgG+IgM H&L preadsorbed	
<b>Target species</b>	Rabbit	
<b>Specificity</b>	Reacts with Rabbit IgG and Rabbit IgM as demonstrated by ELISA; may also react with Ig from other species.	
<b>Tested applications</b>	<b>Suitable for:</b> Flow Cyt, ELISA	
<b>Minimal cross-reactivity</b>	Dog, Goat, Guinea pig, Hamster, Human, Monkey, Mouse, Pig, Rat, Sheep	<a href="#">more details</a>
<b>Immunogen</b>	The details of the immunogen for this antibody are not available.	

### Properties

<b>Form</b>	Liquid	
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.	
<b>Storage buffer</b>	pH: 8.20 Constituent: 100% Borate buffered saline	
	No preservatives or amine-containing buffer salts added.	
<b>Purity</b>	IgG fraction	
<b>Clonality</b>	Monoclonal	
<b>Clone number</b>	SB87a	
<b>Isotype</b>	IgG2a	

### Applications

Our [Abpromise guarantee](#) covers the use of **ab80301** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent dilution.
ELISA		1/4000 - 1/8000.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors