

## Product datasheet

# Recombinant Human Syndecan-1 protein (Tagged) ab83609

### 1 References

#### Overview

---

<b>Product name</b>	Recombinant Human Syndecan-1 protein (Tagged)
<b>Protein length</b>	Protein fragment

#### Description

---

<b>Nature</b>	Recombinant
<b>Source</b>	Escherichia coli

#### Amino Acid Sequence

<b>Species</b>	Human
<b>Amino acids</b>	18 to 218

#### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab83609** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	ELISA Western blot
---------------------	-----------------------

<b>Form</b>	Liquid
-------------	--------

#### Preparation and Storage

---

<b>Stability and Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Preservative: None Constituents: 20% Glycerol, 50mM Tris acetate, 1mM EDTA, pH 7.5 BSA and Azide free
------------------------------	--

#### General Info

---

<b>Function</b>	Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix.
<b>Sequence similarities</b>	Belongs to the syndecan proteoglycan family.
<b>Cellular localization</b>	Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors