

Product datasheet

Anti-ZNF571 antibody ab85844

1 Image

Overview

Product name	Anti-ZNF571 antibody
Description	Rabbit polyclonal to ZNF571
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 144-193 of Human ZNF571 (CRQGFSYLSCLIQHEENHNIEKCSEVKKHRNTFSKKPSYIQHQRIQTGEK) NP_057620. <div style="text-align: right;"> Run BLAST with Run BLAST with </div>
Positive control	HeLa whole cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 2% Sucrose, PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab85844** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 71 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

Application	Abreviews	Notes
-------------	-----------	-------

Target

Function	May be involved in transcriptional regulation.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 17 C2H2-type zinc fingers. Contains 1 KRAB domain.
Cellular localization	Nucleus.

Images



Anti-ZNF571 antibody (ab85844) at 1 µg/ml +
HeLa whole cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 71 kDa

Observed band size: 71 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors