

Product datasheet

Anti-ALDH7A1 antibody ab92880

[2 Images](#)

Overview

Product name	Anti-ALDH7A1 antibody
Description	Mouse polyclonal to ALDH7A1
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human
Immunogen	Full length protein, corresponding to amino acids 1-511 of Human ALDH7A1 (NP_001173.1).
Positive control	Human liver tissue lysate and ALDH7A1 transfected 293T cell lysate

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Storage buffer	Preservative: None Constituents: 1X PBS, pH 7.2
Purity	Protein A purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab92880** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		

Application notes	WB: 1/500 - 1/1000. Predicted molecular weight: 55 kDa. Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.
--------------------------	--


Target

Images



Anti-ALDH7A1 antibody (ab92880) at 1/500 dilution + Human liver tissue lysate at 50 µg

Predicted band size : 55 kDa



All lanes : Anti-ALDH7A1 antibody (ab92880) at 1/500 dilution

Lane 1 : ALDH7A1 transfected 293T cell lysate

Lane 2 : Non-transfected lysate.

Lysates/proteins at 50 µg per lane.

Predicted band size : 55 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors