

Product datasheet

Goat Anti-Human IgG H&L (DyLight® 550) preadsorbed ab98629

1 Image

Overview

<b>Product name</b>	Goat Anti-Human IgG H&L (DyLight® 550) preadsorbed
<b>Host species</b>	Goat
<b>Target species</b>	Human
<b>Tested applications</b>	<b>Suitable for:</b> IHC-P, ICC/IF, Flow Cyt
<b>Minimal cross-reactivity</b>	Mouse, Rat <a href="#">more details</a>
<b>Conjugation</b>	DyLight® 550. Ex: 562nm, Em: 576nm

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C.
<b>Storage buffer</b>	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS
<b>Purity</b>	Immunogen affinity purified
<b>Purification notes</b>	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 550.
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>General notes</b>	DyLight® 550 replaces DyLight® 549.

Applications


Our [Abpromise guarantee](#) covers the use of **ab98629** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
-------------	-----------	-------

Application	Abreviews	Notes
IHC-P		1/50 - 1/500. DyLight® 550 replaces DyLight® 549. The DyLight® 550 excitation is 562nm and the emission is 576nm. DyLight® 550 provides the same orange-to-red fluorescence and photostability as the DyLight® 549 dye. Please refer to the vial in order to check whether the product you received is DyLight® 549 or DyLight® 550. This information however will have no impact on your experiments. For more information please refer to <a href="#">DyLight® 550</a> .
ICC/IF		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

## Images



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

DyLight®-Goat polyclonal Secondary Antibody to Human IgG - H&L (DyLight® 550), pre-adsorbed(ab98629)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors