

Product datasheet

Mouse monoclonal [SA-DA4] Anti-Human IgM mu chain (HRP) ab99737

Overview

Product name	Mouse monoclonal [SA-DA4] Anti-Human IgM mu chain (HRP)
Host species	Mouse
Target species	Human
Specificity	ab99737 reacts with the heavy chain of Human IgM as demonstrated by ELISA; may also react with IgM from other species.
Tested applications	Suitable for: ELISA, WB, Dot blot, IHC-Fr, ICC
Conjugation	HRP

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: None Constituents: 50% PBS, 50% Glycerol, pH 7.4
Purity	Protein G purified
Clonality	Monoclonal
Clone number	SA-DA4
Isotype	IgG1

Applications

Our [Abpromise guarantee](#) covers the use of **ab99737** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/1000 - 1/2000.
WB		Use at an assay dependent dilution.
Dot blot		Use at an assay dependent dilution.

Application	Abreviews	Notes
IHC-Fr		Use at an assay dependent dilution.
ICC		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors