

Product datasheet

Mouse monoclonal [HP6050] Anti-Human IgG3 hinge heavy chain (FITC) ab99827

Overview

Product name	Mouse monoclonal [HP6050] Anti-Human IgG3 hinge heavy chain (FITC)
Host species	Mouse
Target species	Human
Specificity	Reacts with the hinge region of the heavy chain of Human IgG<>.
Tested applications	Suitable for: Dot blot, ELISA, IHC-Fr, ICC/IF, WB
Conjugation	FITC. Ex: 493nm, Em: 528nm

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: 0.1% Sodium Azide Constituents: PBS
Purity	Protein G purified
Purification notes	Purified from ascites.
Clonality	Monoclonal
Clone number	HP6050
Isotype	IgG1
Light chain type	kappa

Applications

Our [Abpromise guarantee](#) covers the use of **ab99827** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Dot blot		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.

Application	Abreviews	Notes
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		Use 1 µg for 10 ⁶ cells.
WB		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors