

Product datasheet

Rat monoclonal [SB87 α] Anti-Rabbit IgG+IgM H&L (HRP) preadsorbed ab99848

Overview

Product name	Rat monoclonal [SB87a] Anti-Rabbit IgG+IgM H&L (HRP) preadsorbed	
Description	Rat monoclonal [SB87a] Anti-Rabbit IgG+IgM H&L (HRP) preadsorbed	
Host species	Rat	
Target species	Rabbit	
Specificity	ab99848 reacts with Rabbit IgG and Rabbit IgM as demonstrated by ELISA. It has been tested by ELISA to ensure minimal cross reactivity to human, goat, mouse, rat, hamster, pig, guinea pig, dog, sheep, and monkey immunoglobulins. Reactivity to species not listed has not been determined.	
Tested applications	Suitable for: ELISA	
Minimal cross-reactivity	Dog, Goat, Guinea pig, Hamster, Human, Monkey, Mouse, Pig, Rat, Sheep	more details
Conjugation	HRP	

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	Preservative: None Constituents: 50% PBS, 50% Glycerol, pH 7.4
Purity	Protein A purified
Clonality	Monoclonal
Clone number	SB87a
Isotype	IgG2a

Applications

Our [Abpromise guarantee](#) covers the use of **ab99848** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		1/4000 - 1/8000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors