


Product datasheet

Anti-pro Guanylin antibody ab14428

Overview

<b>Product name</b>	Anti-pro Guanylin antibody
<b>Description</b>	Rabbit polyclonal to pro Guanylin
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> ELISA, WB, RIA
<b>Species reactivity</b>	<b>Reacts with:</b> Human <b>Predicted to work with:</b> Guinea pig, Pig 
<b>Immunogen</b>	Synthetic peptide: ESVKKLKDLDQEPQ , corresponding to amino acids 13-25 of Human pro Guanylin. <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
<b>Storage buffer</b>	Preservative: None Constituents: Whole Serum
<b>Purity</b>	Whole antiserum
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab14428** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Predicted molecular weight: 13 kDa.

Application	Abreviews	Notes
RIA		Use at an assay dependent dilution.

---

## Target

**Relevance** Guanylin is an endogenous activator of intestinal guanylate cyclase. It stimulates this enzyme through the same receptor binding region as the heat-stable enterotoxins. Guanylin is highly expressed in ileum and colon, found in plasma and belongs to the guanylin family.

**Cellular localization** Secreted protein

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.com/abpromise> or contact our technical team.

---

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors