

Product datasheet

Normal Guinea Pig Serum (sterile) ab7482

Overview

Product name	Normal Guinea Pig Serum (sterile)
Host species	Guinea pig
Tested applications	Suitable for: Blocking, IHC-P, IHC-Fr
General notes	<p>Total Protein Concentration : normal. Strain: English Short Hair. Sex: Mixed sex</p> <p>Raised in: Guinea pig</p> <p>Serums are obtained from non-hemolyzed blood that is collected from healthy and fasted donors. The donors are not medicated in any way and are maintained on an antibiotic free diet. Sterilization is accomplished by Millipore filtration. These serums provide excellent growth promoting components for tissue culture and microbiological organisms.</p> <p>Purity: Whole antiserum</p> <p>Sterilization: This product was aseptically filtered through a 0.22 micron filter into clean, pre-sterilized containers. The product was tested on trypticase soy agar for 24 hours, 48 hours and 72 hours and was found to be negative for bacteria.</p>

Properties

Form	Liquid
Storage instructions	Shipped on Dry Ice. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab7482 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Blocking		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors