

Product datasheet

KLHL4 overexpression 293T lysate (whole cell) ab107454

Overview

Product name KLHL4 overexpression 293T lysate (whole cell)

Tested applications **Suitable for:** WB

Properties

Mycoplasma free Yes

Form Liquid

Storage instructions Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Storage buffer Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)

Background Similarity: Contains 1 BTB (POZ) domain. Contains 6 Kelch repeats. Tissue specificity: Expressed in adult fibroblasts and in a range of fetal tissues including tongue, palate, and mandible.

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab107454 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors