

Product datasheet

PL6 overexpression 293T lysate (whole cell) ab94331

2 Images

Overview

Product name	PL6 overexpression 293T lysate (whole cell)
General notes	ab94331 is a 293T cell transfected lysate in which Human PL6 has been transiently over-expressed using a pCMV-PL6 plasmid. The lysate is provided in 1X Sample Buffer. <b>Note:</b> For more details about how the transfected lysate was prepared <a href="#">view preparation notes</a>
Tested applications	<b>Suitable for:</b> WB

Properties


Mycoplasma free	Yes
Form	Liquid
Storage instructions	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Storage buffer	pH: 6.80 Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol
Background	PL6 is an integral membrane protein with signal transduction activity.

Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab94331 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		Use at an assay dependent dilution.

Images



SDS-PAGE - PL6 overexpression 293T lysate  
(whole cell) (ab94331)

ab94331 at 15µg/lane on an SDS-PAGE gel.



Western blot - PL6 overexpression 293T lysate  
(whole cell) (ab94331)

**All lanes :** Anti-PL6 antibody (**ab56868**) at 1/500 dilution

**Lane 1 :** PL6 overexpression 293T lysate (whole cell) (ab94331)

**Lane 2 :** 293T Non Transfected Lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

**All lanes :** Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors