


Product datasheet

Anti-alpha 1 Spectrin antibody [EPR9300] ab139403

Recombinant RabMAb

[3 References](#) [3 Images](#)

Overview

Product name	Anti-alpha 1 Spectrin antibody [EPR9300]
Description	Rabbit monoclonal [EPR9300] to alpha 1 Spectrin
Host species	Rabbit
Tested applications	Suitable for: Flow Cyt (Intra), WB Unsuitable for: IHC-P or IP
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat 
Immunogen	Synthetic peptide within Human alpha 1 Spectrin aa 2300-2400. The exact sequence is proprietary.
Positive control	K562, Human fetal heart and fetal liver lysates; K562 cells
General notes	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at -20°C.
Storage buffer	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 9% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA, 50% Tissue culture supernatant
Purity	Tissue culture supernatant
Clonality	Monoclonal

Clone number	EPR9300
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab139403 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
Flow Cyt (Intra)		1/100 - 1/500. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
WB		1/1000 - 1/10000. Predicted molecular weight: 280 kDa.

Application notes Is unsuitable for IHC-P or IP.

Target

Function	Spectrin is the major constituent of the cytoskeletal network underlying the erythrocyte plasma membrane. It associates with band 4.1 and actin to form the cytoskeletal superstructure of the erythrocyte plasma membrane.
Involvement in disease	<p>Defects in SPTA1 are the cause of elliptocytosis type 2 (EL2) [MIM:130600]. EL2 is a Rhesus-unlinked form of hereditary elliptocytosis, a genetically heterogeneous, autosomal dominant hematologic disorder. It is characterized by variable hemolytic anemia and elliptical or oval red cell shape.</p> <p>Defects in SPTA1 are a cause of hereditary pyropoikilocytosis (HPP) [MIM:266140]. HPP is an autosomal recessive disorder characterized by hemolytic anemia, microspherocytosis, poikilocytosis, and an unusual thermal sensitivity of red cells.</p> <p>Defects in SPTA1 are the cause of spherocytosis type 3 (SPH3) [MIM:270970]; also known as hereditary spherocytosis type 3 (HS3). Spherocytosis is a hematologic disorder leading to chronic hemolytic anemia and characterized by numerous abnormally shaped erythrocytes which are generally spheroidal. SPH3 is characterized by severe hemolytic anemia. Inheritance is autosomal recessive.</p>
Sequence similarities	<p>Belongs to the spectrin family.</p> <p>Contains 3 EF-hand domains.</p> <p>Contains 1 SH3 domain.</p> <p>Contains 21 spectrin repeats.</p>
Cellular localization	Cytoplasm > cytoskeleton. Cytoplasm > cell cortex.

Images



Western blot - Anti-alpha 1 Spectrin antibody
[EPR9300] (ab139403)

All lanes : Anti-alpha 1 Spectrin antibody [EPR9300] (ab139403)
at 1/1000 dilution

Lane 1 : K562 lysate

Lane 2 : Human fetal heart lysate


Lane 3 : Human fetal liver lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : HRP labelled goat anti-rabbit at 1/2000 dilution

Predicted band size: 280 kDa



Flow Cytometry (Intracellular) - Anti-alpha 1 Spectrin
antibody [EPR9300] (ab139403)

Intracellular flow cytometric analysis of permeabilized K562 cells
labelling alpha 1 Spectrin with ab139403 at 1/100 dilution (red) or a
rabbit IgG (negative) (green).

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-alpha 1 Spectrin antibody [EPR9300]
(ab139403)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors