

Product datasheet

Anti-CTRP6 antibody ab36898

3 Images

Overview

Product name	Anti-CTRP6 antibody
Description	Rabbit polyclonal to CTRP6
Host species	Rabbit
Tested applications	Suitable for: ICC/IF, WB
Species reactivity	Reacts with: Human
Immunogen	A 15 amino acid peptide from near the carboxy-terminus of human CTRP6.
Positive control	HeLa whole cell lysate (ab150035)
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications


The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab36898 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		Use a concentration of 10 µg/ml.
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 29 kDa.

Target

Sequence similarities	Contains 1 C1q domain. Contains 1 collagen-like domain.
Cellular localization	Secreted.

Images



Lane 1 : Anti-CTRP6 antibody (ab36898) at 1 µg/ml

Lane 2 : Anti-CTRP6 antibody (ab36898) at 2 µg/ml

All lanes : Hela cell lysate.


Predicted band size: 29 kDa

Observed band size: 42 kDa

These proteins often migrate in SDS-PAGE at positions other than their predicted size.



ab36898 at 10µg/ml staining CTRP6 in Hela cells by ICC/IF



Immunofluorescence of CTRP6 in Hela cells using ab36898 at 20 ug/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors