

Product datasheet

FITC Mouse IgG2α, monoclonal - Isotype control ab81197

4 References

Overview

Product name	FITC Mouse IgG2α, monoclonal - Isotype control
Conjugation	FITC. Ex: 493nm, Em: 528nm
Conjugation notes	Excitation: 488 nm. Peak emission: 525 nm.
Tested applications	Suitable for: Flow Cyt
General notes	<p>Reagents should be brought to room temperature ($22 \pm 3^{\circ}\text{C}$) before use. Protect cells from light sources during incubation with antibodies and prior to analysis.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C . Store at $+4^{\circ}\text{C}$ short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. Store In the Dark.
Storage buffer	<p>pH: 7.3</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituents: 0.4% BSA, 1.74% Sodium chloride, 0.328% Sodium phosphate</p>
Purity	IgG fraction
Clonality	Monoclonal
Clone number	TIB13
Isotype	IgG2α

Applications

The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab81197 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Flow Cyt		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors