

# Anti-Ganglioside OAcGD3 antibody [7H2] ab6218

### Overview

<b>Product name</b>	Anti-Ganglioside OAcGD3 antibody [7H2]
<b>Description</b>	Mouse monoclonal [7H2] to Ganglioside OAcGD3
<b>Host species</b>	Mouse
<b>Specificity</b>	This antibody reacts specifically with O-AcGD3 ganglioside: OAcNeuAcá2-8NeuAcá2-3Galâ1-4Galcâ1-1Cer.
<b>Tested applications</b>	<b>Suitable for:</b> WB, ELISA, Thin Layer Chromatography
<b>General notes</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### Properties

<b>Form</b>	Lyophilized:Reconstitute with 1 ml sterile water.
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
<b>Storage buffer</b>	Constituent: Tissue culture supernatant
<b>Purity</b>	Tissue culture supernatant
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	7H2
<b>Isotype</b>	IgG3
<b>Light chain type</b>	kappa

### Applications

**The Abpromise guarantee** Our **Abpromise guarantee** covers the use of ab6218 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/10 - 1/100.
ELISA		1/10 - 1/100.
Thin Layer Chromatography		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors