

Product datasheet

Anti-hapten 4-hydroxy-3-nitrophenyl acetyl (NP) antibody [B1-8] ab206539

Recombinant

[1 References](#) [1 Image](#)

Overview

| | |
|----------------------------|---|
| Product name | Anti-hapten 4-hydroxy-3-nitrophenyl acetyl (NP) antibody [B1-8] |
| Description | Rat monoclonal [B1-8] to hapten 4-hydroxy-3-nitrophenyl acetyl (NP) |
| Host species | Rat |
| Tested applications | Suitable for: WB, Flow Cyt |
| Species reactivity | Reacts with: Species independent |
| Immunogen | Recombinant fragment. Hapten-protein conjugated NP12-chicken gammaglobulin (CG) |

Properties

| | |
|-----------------------------|---|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. |
| Storage buffer | Preservative: 0.02% Proclin 300 Constituent: 99% PBS |
| Purity | Protein A purified |
| Clonality | Monoclonal |
| Clone number | B1-8 |
| Isotype | IgG2c |
| Light chain type | lambda |

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab206539 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|--|
| WB | | Use at an assay dependent concentration. |
| Flow Cyt | | Use at an assay dependent concentration. |

Images

Why choose a recombinant antibody?




Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-hapten 4-hydroxy-3-nitrophenyl acetyl (NP) antibody [B1-8] (ab206539)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors