


Product datasheet

Anti-Jagged1 antibody [EPR4290] ab109536

KO VALIDATED Recombinant RabMAB

★★★★★ [4 Abreviews](#) [52 References](#) [5 Images](#)

Overview

Product name	Anti-Jagged1 antibody [EPR4290]
Description	Rabbit monoclonal [EPR4290] to Jagged1
Host species	Rabbit
Tested applications	Suitable for: WB Unsuitable for: ICC/IF or IHC-P
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat 
Immunogen	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
Positive control	WB: HepG2 and HUVEC whole cell lysate, NIH:OVCA9-3 cell lysate.
General notes	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production <p>For more information see here.</p> <p>Our RabMAB[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAB[®] patents.</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at -20°C.
Storage buffer	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 40% Glycerol, 59% PBS, 0.05% BSA
Purity	Protein A purified
Clonality	Monoclonal
Clone number	EPR4290
Isotype	IgG

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab109536 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB	★★★★★ (2)	1/1000 - 1/10000. Predicted molecular weight: 134 kDa.

Application notes Is unsuitable for ICC/IF or IHC-P.

Target

Function Ligand for multiple Notch receptors and involved in the mediation of Notch signaling. May be involved in cell-fate decisions during hematopoiesis. Seems to be involved in early and late stages of mammalian cardiovascular development. Inhibits myoblast differentiation (By similarity). Enhances fibroblast growth factor-induced angiogenesis (in vitro).

Tissue specificity Widely expressed in adult and fetal tissues. In cervix epithelium expressed in undifferentiated subcolumnar reserve cells and squamous metaplasia. Expression is up-regulated in cervical squamous cell carcinoma. Expressed in bone marrow cell line HS-27a which supports the long-term maintenance of immature progenitor cells.

Involvement in disease Defects in JAG1 are the cause of Alagille syndrome type 1 (ALGS1) [MIM:118450]. Alagille syndrome is an autosomal dominant multisystem disorder defined clinically by hepatic bile duct paucity and cholestasis in association with cardiac, skeletal, and ophthalmologic manifestations. There are characteristic facial features and less frequent clinical involvement of the renal and vascular systems.


Defects in JAG1 are a cause of tetralogy of Fallot (TOF) [MIM:187500]. TOF is a congenital heart anomaly which consists of pulmonary stenosis, ventricular septal defect, dextroposition of the aorta (aorta is on the right side instead of the left) and hypertrophy of the right ventricle. This condition results in a blue baby at birth due to inadequate oxygenation. Surgical correction is emergent.

Sequence similarities Contains 1 DSL domain.
Contains 15 EGF-like domains.

Developmental stage Expressed in 32-52 days embryos in the distal cardiac outflow tract and pulmonary artery, major arteries, portal vein, optic vesicle, otocyst, branchial arches, metanephros, pancreas, mesocardium, around the major bronchial branches, and in the neural tube.

Cellular localization Membrane.

Images



All lanes : Anti-Jagged1 antibody [EPR4290] (ab109536) at 1/10000 dilution

Lane 1 : HepG2 (human hepatocellular carcinoma epithelial cell), whole cell lysate

Lane 2 : HUVEC (human umbilical vein endothelial cell), whole cell lysate

Lane 3 : Jurkat (human T cell leukemia T lymphocyte), whole cell lysate

Lane 4 : MOLT-4 (human lymphoblastic leukemia T lymphoblast), whole cell lysate

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

Predicted band size: 134 kDa

Observed band size: 200, 23 kDa

Exposure time: 180 seconds


Blocking buffer and concentration: 5% NFDm/TBST

Diluting buffer and concentration: 5% NFDm /TBST

Negative controls: Jurkat and MOLT-4 whole cell lysate (PMID: 30231940)

200-kDa full length and 23-kDa C-terminal JAG1 are observed. The molecular weights are consistent with what has been described in the literature (PMID: 30890522, 30890522).

This blot was developed using a high sensitivity ECL substrate.



Western blot - Anti-Jagged1 antibody [EPR4290] (ab109536)

All lanes : Anti-Jagged1 antibody [EPR4290] (ab109536) at 1/10000 dilution

Lane 1 : Wild-type HAP1 whole cell lysate

Lane 2 : JAG1 knockout HAP1 whole cell lysate


Lysates/proteins at 20 µg per lane.

Predicted band size: 134 kDa

Observed band size: 200 kDa

Lanes 1 - 2: Merged signal (red and green). Green - ab109536 observed at 200 kDa. Red - loading control, **ab9484**, observed at 37 kDa.

ab109536 was shown to recognize in wild-type HAP1 cells as signal was lost at the expected MW in JAG1 knockout cells. Additional cross-reactive bands were observed in the wild-type and knockout cells. Wild-type and JAG1 knockout samples were subjected to SDS-PAGE. Ab109536 and **ab9484** (Mouse anti GAPDH loading control) were incubated overnight at 4°C at 1/10000 dilution and 1/20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preabsorbed **ab216773** and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed **ab216776** secondary antibodies at 1/20000 dilution for 1 hour at room temperature before imaging.



Western blot - Anti-Jagged1 antibody [EPR4290] (ab109536)

Anti-Jagged1 antibody [EPR4290] (ab109536) at 1/1000 dilution (Purified) + HepG2 at 10 µg

Secondary


HRP goat anti-rabbit (H+L) at 1/1000 dilution

Predicted band size: 134 kDa

Observed band size: 200 kDa

Blocking buffer: 5% NFDM/TBST

Dilution buffer: 5% NFDM/TBST



Western blot - Anti-Jagged1 antibody [EPR4290]
(ab109536)

All lanes : Anti-Jagged1 antibody [EPR4290] (ab109536) at 1/10000 dilution (Unpurified)

Lane 1 : HepG2 cell lysate

Lane 2 : NIH:OVCA3 cell lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : HRP-labelled goat anti-rabbit at 1/2000 dilution

Predicted band size: 134 kDa

Observed band size: 200 kDa

Why choose a recombinant antibody?

<p>Research with confidence Consistent and reproducible results</p>	<p>Long-term and scalable supply Recombinant technology</p>
<p>Success from the first experiment Confirmed specificity</p>	<p>Ethical standards compliant Animal-free production</p>

Anti-Jagged1 antibody [EPR4290] (ab109536)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors