

# Anti-Mumps nucleoprotein antibody [5E3] ab9879

[2 References](#) [1 Image](#)

### Overview

<b>Product name</b>	Anti-Mumps nucleoprotein antibody [5E3]
<b>Description</b>	Mouse monoclonal [5E3] to Mumps nucleoprotein
<b>Host species</b>	Mouse
<b>Tested applications</b>	<b>Suitable for:</b> Indirect ELISA, WB
<b>Species reactivity</b>	<b>Reacts with:</b> Mumps virus
<b>Immunogen</b>	Recombinant mumps NP. Gloucester strain (Glouc 1/UK96)
<b>General notes</b>	<p>This product was changed from ascites to tissue culture supernatant on 28/11/2017. Please note that the dilutions may need to be adjusted accordingly.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.2 Preservative: 0.1% Sodium azide
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	5E3
<b>Myeloma</b>	NS0/1
<b>Isotype</b>	IgG2a
<b>Light chain type</b>	unknown

## Applications

### The Abpromise guarantee


Our **Abpromise guarantee** covers the use of ab9879 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
Indirect ELISA		1/1000 - 1/5000.
WB		1/500 - 1/1000.

## Target

## Images



Western blot of yeast-derived mumps NP using ab9879 at a concentration of 0.5 µg/ml.

Western blot of yeast-derived mumps NP using ab9879 at a concentration of 0.5 µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors