

Product datasheet

Anti-prekallikrein + Kallikrein immobilized antibody [10G3-42] ab185493

Overview

Product name	Anti-prekallikrein + Kallikrein immobilized antibody [10G3-42]
Description	Mouse monoclonal [10G3-42] to prekallikrein + Kallikrein immobilized
Host species	Mouse
Immunogen	Full length native protein (purified) corresponding to Human prekallikrein + Kallikrein immobilized aa 1 to the C-terminus. Purified from Hybridoma cell culture. Database link: P06870

 [Run BLAST with](#)

 [Run BLAST with](#)

General notes

ab185493 is an immobilized monoclonal antibody to Human prekallikrein and kallikrein.

Supplied as a resin; 1/1 slurry.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

Properties

Form	Solid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.50 Preservative: 0.02% Sodium azide Constituents: 0.58% Sodium chloride, 1.21% Tris
Purity	Protein G purified
Clonality	Monoclonal
Clone number	10G3-42
Isotype	IgG1
Light chain type	kappa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors