

Product datasheet

Anti-Syndecan-1 antibody [B-A38], prediluted ab1410

2 References

Overview

Product name	Anti-Syndecan-1 antibody [B-A38], prediluted
Description	Mouse monoclonal [B-A38] to Syndecan-1, prediluted
Host species	Mouse
Tested applications	Suitable for: IHC-Fr, IHC-P
Species reactivity	Reacts with: Human
Immunogen	Tissue, cells or virus corresponding to Syndecan-1. U266 cell line derived from Hodgkins disease
Positive control	Tonsil
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.30 Constituents: PBS, Tissue culture supernatant
Clonality	Monoclonal
Clone number	B-A38
Myeloma	unknown
Isotype	IgG1
Light chain type	unknown

Applications

Applications

The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab1410 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-Fr		
IHC-P		

Application notes

IHC:

This antibody is pre-diluted for use with an ABC system.

For formalin fixed paraffin embedded sections, perform heat mediated antigen retrieval before commencing with IHC.

Target

Function

Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix.

Sequence similarities

Belongs to the syndecan proteoglycan family.

Cellular localization

Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors