

Anti-TGE Virus peplomer antibody [1Q17] ab2599

Overview

Product name	Anti-TGE Virus peplomer antibody [1Q17]
Description	Mouse monoclonal [1Q17] to TGE Virus peplomer
Host species	Mouse
Specificity	This antibody is specific for TGE (Transmissible gastroenteritis) virus peplomer antigen.
Tested applications	Suitable for: WB, ELISA
Species reactivity	Reacts with: Transmissible gastroenteritis
Immunogen	Tissue, cells or virus corresponding to TGE Virus peplomer. Purified TGE virus
Epitope	The exact epitope recognised by this antibody has not been mapped.
General notes	<p>The peplomer protein mediates the binding of virions to the host cell receptor and is also involved with membrane fusion.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.20 Preservative: 0.01% Sodium azide
Primary antibody notes	The peplomer protein mediates the binding of virions to the host cell receptor and is also involved with membrane fusion.
Clonality	Monoclonal
Clone number	1Q17
Myeloma	unknown

Isotype	IgG1
Light chain type	unknown

Applications

The Abpromise guarantee Our **Abpromise guarantee** covers the use of ab2599 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000.
ELISA		1/10000.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors