

Product datasheet

Anti-ZNF704 antibody ab121804

[1 References](#) [2 Images](#)

Overview

Product name	Anti-ZNF704 antibody
Description	Rabbit polyclonal to ZNF704
Host species	Rabbit
Tested applications	Suitable for: WB, ICC/IF
Species reactivity	Reacts with: Human
Immunogen	antigen: TIPGSAKFTP NGSSFISWQ SPPVTFTGIP VSPTHHPVGT GEQRQHAHTV LSSPPRGTVS LRKPRGEGKK CRKVYGMENR DMWCT, corresponding to amino acids 315-399 of Human ZNF704 (Q6ZNC4).

[Run BLAST with Expasy](#) [Run BLAST with NCBI](#)

Positive control Human small intestine tissue

General notes The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
Storage buffer	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee

Our **Abpromise guarantee** covers the use of ab121804 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		Use a concentration of 0.04 - 0.4 µg/ml.
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment

Target

Sequence similarities


Contains 1 C2H2-type zinc finger.

Images



Western blot - Anti-ZNF704 antibody (ab121804)

Anti-ZNF704 antibody (ab121804) at 0.04 µg/ml (Anti-ZNF704 primary antibody on human cell-line SH-SY5Y)



Immunocytochemistry/ Immunofluorescence - Anti-ZNF704 antibody (ab121804)

Immunofluorescent staining of Human cell line A-431 shows positivity in nucleus but not nucleoli. Recommended concentration of ab121804 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors