

Product datasheet

Human Thymulin peptide ab42294

Description

Product name	Human Thymulin peptide
Purity	> 95 % SDS-PAGE. Purified by HPLC.
Animal free	No
Nature	Synthetic
Species	Human
Sequence	EAKSQGGSN

Specifications

Our **Abpromise guarantee** covers the use of **ab42294** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Blocking
Form	Lyophilized

Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Reconstitution	Resuspend in 0.1% acetic acid for required concentration

General Info

Relevance	Thymulin, a nonapeptide secreted by the thymus, is an essential hormone for T lymphocyte differentiation and function.
Cellular localization	Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors