

Product datasheet

Recombinant Hepatitis C Virus NS4 protein (Rhodamine) **ab49031**

Description

Product name	Recombinant Hepatitis C Virus NS4 protein (Rhodamine)
Purity	> 95 % SDS-PAGE.
Expression system	Escherichia coli
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Amino acids	1658 to 1863
Conjugation	Rhodamine. Ex: 550nm, Em: 570nm

Specifications

Our **Abpromise guarantee** covers the use of **ab49031** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Form	Liquid
-------------	--------

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle. pH: 7.2 Constituents: 48% Urea, 0.078% Beta mercaptoethanol, 0.316% Tris HCl
------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

General Info

Relevance	NS4 A/B are two of the seven nonstructural (NS) proteins making up the HCV polyprotein. The NS proteins recruit the viral genome into an RNA replication complex, which is associated with rearranged cytoplasmic membranes. NS4A acts as a cofactor with the NS3 serine protease and stabilizes its folding. The NS3-NS4A complex is essential for the activation of the latter and allows membrane anchorage of NS3. HCV is a positive, single-stranded RNA virus in the Flaviviridae family. The genome is approximately 10,000 nucleotides and encodes a single polyprotein of
------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

about 3,000 amino acids. HCV is responsible for a large proportion of worldwide chronic viral hepatitis. Most of these infections develop into chronic hepatitis, which often progresses to liver cirrhosis and hepatocellular carcinoma. At present, (unlike hepatitis A and B), there is no vaccine to prevent hepatitis C infection. The hepatitis C virus (HCV) nonstructural protein 4B (NS4B) is a relatively hydrophobic 27-kDa protein. The 4A protein has a molecular weight of 6 kDa.

Cellular localization

Endoplasmic reticulum; multi-pass membrane protein.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors