

Recombinant Human Acyloxyacyl Hydrolase protein

ab157892

1 Image

Description	
Product name	Recombinant Human Acyloxyacyl Hydrolase protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<div>MQSPWKILTVAPLFLLLSLQSSASPANDDQSRPSLSNGHT</div> <div>CVGCVLVVSV</div> <div>IEQLAQVHNSTVQASMERLCSYLPEKLFLKTTCYLVIDKFG</div> <div>SDIILKLLSA</div> <div>DMNADVVCHTLEFCKQNTGQPLCHLYPLPKETWKFTLQK</div> <div>ARQIVKKSPIL</div> <div>KYSRSGSDICSLPVLAKICQIKLAMEQSVPFKDVDSDKY</div> <div>SVFPTLRGYH</div> <div>WRGRDCNDSDSVYPGRRPNNWDVHQDSNCNGWGV</div> <div>PKDGVPEKKFCE</div> <div>GSQPRGIILLGDSAGAHFHISPEWITASQMSLNSFINLPTAL</div> <div>TNELDWPQ</div> <div>LSGATGFLDSTVGIKEKSYLRLWKRNHCHNRDYQNISRNG</div> <div>ASSRNLKKF</div> <div>IESLSRNKVLDPAMVYAMIGNDVCSGKSDPVPAMTTPEK</div> <div>LYSNVMQTL</div> <div>KHLNSHLPNGSHVILYGLPDGTFLWDNLHNRYHPLGQLNK</div> <div>DMTYAQLYSF</div> <div>LNCLQVSPCHGWMSSNKTTLTTERAEQLSNTLKKAAS</div> <div>EKFTNFNLFY</div> <div>MDFAFHEIIQEWQKRGGQPWQLIEPVDGFHPNEVALLLA</div> <div>DHFWKKVQLQ WPQILGKENPFNPQIKQVFGDQGGH</div>
Amino acids	1 to 575
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab157892** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA
	Western blot
Form	Liquid
Additional notes	


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
	pH: 8.00
	Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function	Removes the secondary (acyloxyacyl-linked) fatty acyl chains from the lipid A region of bacterial lipopolysaccharides.
Sequence similarities	Contains 1 saposin B-type domain.
Post-translational modifications	Both subunits contain a number of cysteine residues that may form disulfide bridges.
Cellular localization	Secreted.

Images



ab157892 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors