

Product datasheet

Recombinant Human ADAL protein (denatured)
ab180314

1 Image

Description	
Product name	Recombinant Human ADAL protein (denatured)
Purity	> 80 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>Q6DHV7-2</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSHMGSEFMIEAEEQQPCKTD FYSELPKVELHA HLNGSISSTMTMKLIAQKPDLDQMTVIDKGKKRTLEEC FQMFQTIHQ LTSSPEDILMVTKDVIKEFADDGVKYLELRSTPRRENATG MTKKTYVESI LEGIKQSKQENLDIDVRYLIAVDRRGGPLVAKETVKLAEEF FLSTEGTVL GLDLSGDPTVGQAKDFLEPLLEAKKAGLKLALHLSEIPNQ KKETQILLDL LPDRIGHGTFLNSGEGGSLDLVDFVRQHRIPLGKAWSFRS SR
Predicted molecular weight	33 kDa including tags
Amino acids	1 to 267
Tags	His tag N-Terminus
Additional sequence information	NP_001012987
Description	Recombinant Human ADAL protein

Specifications

Our Abpromise guarantee covers the use of **ab180314** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00


Constituents: 0.32% Tris HCl, 2.4% Urea, 10% Glycerol (glycerin, glycerine)

General Info

Function Putative nucleoside deaminase. May catalyze the hydrolytic deamination of adenosine or some similar substrate and play a role in purine metabolism.

Sequence similarities Belongs to the adenosine and AMP deaminases family.

Images



3ug by SDS-PAGE under reducing conditions and visualized by coomassie blue stain.

SDS-PAGE - Recombinant Human ADAL protein
(denatured) (ab180314)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors