

Recombinant human ANGPTL7 protein ab109158

Description

Product name	Recombinant human ANGPTL7 protein
Biological activity	Activates ERK phosphorylation in THP1 cells.
Purity	> 90 % SDS-PAGE. ab109158 is an 0.2 µm filtered solution.
Endotoxin level	< 0.100 Eu/µg
Expression system	HEK 293 cells
Accession	<u>O43827</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Predicted molecular weight	45 kDa including tags
Amino acids	1 to 346
Tags	DDDDK tag C-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab109158** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituent: PBS This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------------------------	---

General Info

Tissue specificity	Highly and specifically expressed in the cornea where is confined to the stromal layer.
Sequence similarities	Contains 1 fibrinogen C-terminal domain.
Cellular localization	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors