

Recombinant Human BET3 protein ab124564

1 Image

Description	
Product name	Recombinant Human BET3 protein
Purity	> 90 % SDS-PAGE. ab124564 was purified by conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>O43617</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MSRQANRGTE SKKMSSELF LTYGALVTQL CKDYENDEDV NKQLDKMGFN IGVRLIEDFL ARSNVGRCHD FRETADVIAK VAFKMYLGIT PSITNWSPAG DEFSLILENN PLVDFVELPD NHSSLYSNL LCGVLRGALE MVQMAVEAKF VQDTLKGDGV TEIRMRFIRR IEDNLPAGEE
Predicted molecular weight	22 kDa including tags
Amino acids	1 to 180
Tags	His tag N-Terminus

Specifications	
Our Abpromise guarantee covers the use of ab124564 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF
Form	Liquid
Additional notes	Previously labelled as TRAPPC3.

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

General Info

Function

May play a role in vesicular transport from endoplasmic reticulum to Golgi.


Sequence similarities

Belongs to the TRAPP small subunits family. BET3 subfamily.

Cellular localization

Golgi apparatus > cis-Golgi network. Endoplasmic reticulum.

Images



ab124564 (3µg) visualised by SDS-PAGE (15%).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors