

## Product datasheet

### Recombinant Human CCN3 protein ab79564

1 Image

#### Description

<b>Product name</b>	Recombinant Human CCN3 protein
<b>Purity</b>	> 90 % SDS-PAGE.
<b>Expression system</b>	Escherichia coli
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Human
<b>Sequence</b>	KGKKCLRTKKSLKAIHLQFKNCTSLHTYKPRFCGVCSDGR CCTPHNTKI QAEFQCSPGQVKKPVMVIGTCTCHTNCPKNNEAFLQELE LKTTRGKM
<b>Amino acids</b>	260 to 357

#### Specifications

Our **Abpromise guarantee** covers the use of **ab79564** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	Western blot
	ELISA
	SDS-PAGE
<b>Form</b>	Liquid


#### Preparation and Storage

<b>Stability and Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituent: 2.68% PBS
------------------------------	---

#### General Info

<b>Function</b>	Immediate-early protein likely to play a role in cell growth regulation.
<b>Tissue specificity</b>	Expressed in bone marrow, thymic cells and nephroblastoma. Increased expression in Wilms tumor of the stromal type.
<b>Sequence similarities</b>	<p>Belongs to the CCN family.</p> <p>Contains 1 CTCK (C-terminal cystine knot-like) domain.</p> <p>Contains 1 IGFBP N-terminal domain.</p> <p>Contains 1 TSP type-1 domain.</p> <p>Contains 1 VWFC domain.</p>
<b>Cellular localization</b>	Secreted.

## Images



ab79564 on 15% SDS PAGE gel, Coomassie blue staining.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors