

Recombinant Human CD3 epsilon protein ab114153

1 Image

Description

Product name	Recombinant Human CD3 epsilon protein
Expression system	Wheat germ
Accession	<u>P07766</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	DGNEEMGGITQTPYKVSISGTTVILTCPQYPGSEILWQRND KNIGGEDDD KNIGSDEDHLSLKEFSELEQSGYYVCYPRGSKPEDANFYL YLRARVCENC MEMDVMSVATIVIDICITGGLLLL VYWSKNRKAKAKPVT RGAGAGGRQ RGQNKERPPPVNPDPYEPIRKGQRDLYSGLNQRR
Predicted molecular weight	46 kDa including tags
Amino acids	23 to 207
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab114153** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot SDS-PAGE
Form	Liquid

Preparation and Storage


Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
------------------------------	---

Constituents: 0.3% Glutathione, 0.79% Tris HCl

General Info

Function	The CD3 complex mediates signal transduction.
Sequence similarities	Contains 1 Ig-like (immunoglobulin-like) domain. Contains 1 ITAM domain.
Cellular localization	Membrane.

Images



12.5% SDS-PAGE analysis of ab114153, stained with Coomassie Blue.

SDS-PAGE - Recombinant Human CD3 epsilon protein (ab114153)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors