

Product datasheet

Recombinant Human CNOT8 protein ab167839

1 Image

Description

Product name	Recombinant Human CNOT8 protein	
Purity	> 85 % SDS-PAGE. ab167839 was purified using conventional chromatography techniques.	
Expression system	Escherichia coli	
Accession	<u>Q9UFF9</u>	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Human	
Sequence	MGSSHHHHHHSSGLVPRGSHMGSMMPAALVENSQVICEV WASNLEEEMRKI REIVLSYSYIAMDTEFPGVVVRPIGEFRSSIDYQYQLLR CNV DLLKIIQL GLTFTNEKGEYPSGINTWQFNFKFNLTEDMYSQDSIDLLAN SGLQFQKHE EEGIDTLHFAELLMTSGVVLC DNVKWLSFHSGYDFGYMV KLLTDSRLPEE EHEFFHILNFFPSIYDVKYLKMKCKNLKGG LQE VADQLDL QRIGRQHQA GSDSLLTGMAFFRMKELFFEDSIDD AKYCGRLYGLGTGVA QKQNE DVDSA QEKMSILAIINMQQ	
Predicted molecular weight	36 kDa including tags	
Amino acids	1 to 292	
Tags	His tag N-Terminus	

Specifications

Our **Abpromise guarantee** covers the use of **ab167839** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
	Mass Spectrometry
Mass spectrometry	MALDI-TOF

Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00
Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.88% Sodium chloride


General Info

Function Ubiquitous transcription factor required for a diverse set of processes. The CCR4-NOT complex functions as general transcription regulation complex.

Sequence similarities Belongs to the CAF1 family.

Cellular localization Cytoplasm. Nucleus.

Images



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors