

Product datasheet

Recombinant Human CutA protein ab95345

1 Image

Description

Product name	Recombinant Human CutA protein
Purity	> 95 % SDS-PAGE. ab95345 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p> MRLLLLPRVL LTMASGSPPT QPSPASDSGS GYVPGSVSAA FVTC PNEKVA KEIARAVVEK RLAACVNLI P QITSIYEWKG KIEEDSEVLM MIKTQSSLVP ALTDFVRSVH PYEVAEVI AL PVEQGNFPYL QWVRQVTE SV SDSITVLP LE HHHHHH </p>

Specifications

Our **Abpromise guarantee** covers the use of **ab95345** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage


Stability and Storage	<p>Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.</p> <p>pH: 8.00</p> <p>Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)</p>
------------------------------	--

General Info

Function	May form part of a complex of membrane proteins attached to acetylcholinesterase (AChE).
-----------------	--

Tissue specificity	Ubiquitous. Widely expressed in brain.
Sequence similarities	Belongs to the CutA family.
Post-translational modifications	O-glycosylated.

Images



15% SDS-PAGE analysis of 3µg ab95345.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors