

Recombinant Human DCUN1D2 protein ab113157

1 Image

Description	
Product name	Recombinant Human DCUN1D2 protein
Purity	> 95 % SDS-PAGE. ab113157 was purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>Q6PH85</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSH MHKLKSSQKDKVRQFMA CTQAGERTAIYCL TQNEWRLDEATDSFFQNPDSLHRESMRNAVDKKKLERLY GRYKDPQDENK IGVDGIQQFCDDLSLDPASISVLVIAWKFRAATQCEFSRKE FLDGMTELG CDSMEKLKALLPRLEQELKDTAKFKDFYQFTTFKKNPG QKGLDLEMAVA YWKLVLSGRFKFLDLWNTFLMEHHKRSIPRDTWNLLDFFG NMIADDMSNY DEEGAWPVLIDDFVEYARPVVTGGKRSLF
Predicted molecular weight	32 kDa including tags
Amino acids	1 to 259
Tags	His tag N-Terminus

Specifications	
Our Abpromise guarantee covers the use of ab113157 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 30% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

General Info

Function


May contribute to neddylation of cullin components of SCF-type E3 ubiquitin ligase complexes. Neddylation of cullins play an essential role in the regulation of SCF-type complexes activity.

Sequence similarities

Contains 1 DCUN1 domain.

Contains 1 UBA-like domain.

Images



15% SDS PAGE analysis of 3 µg ab113157.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors