

Product datasheet

Recombinant Human DPL protein ab139778

1 Image

Description

Product name	Recombinant Human DPL protein
Purity	> 95 % SDS-PAGE. ab139778 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>Q9UKY0</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSRGIKHRI KWNRKALPST AQITEAQVAE NRPGAFIKQG RKLDIDFGAE GNRYEANYW QFPDGIHYNG CSEANVTKEA FVTGCINATQ AANQGEFQKP DNKLHQQLW RLVQELCSLK HCEFWLERG
Predicted molecular weight	17 kDa including tags
Amino acids	27 to 152
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab139778** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid
Additional notes	This product was previously labelled as Doppel

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.50

Constituents: 0.02% DTT, 0.32% Tris HCl, 30% Glycerol (glycerin, glycerine), 1.17% Sodium chloride

General Info

Tissue specificity

Expressed in testis.


Sequence similarities

Belongs to the prion family.

Cellular localization

Cell membrane.

Images



15% SDS-PAGE analysis of ab139778 (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors