

Recombinant Human FAM19A2 protein ab54415

Description

Product name	Recombinant Human FAM19A2 protein
Purity	> 95 % SDS-PAGE. Purity : Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin Level : < 0.1 ng per µg (1EU/µg).
Endotoxin level	< 0.100 Eu/µg
Expression system	Escherichia coli
Accession	<u>Q8N3H0</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	ANHHKAHHVKTGTCEVVALHRCCNKNKIEERSQTVKCSC FPGQVAGTTRA APSCVDASIVEQKWWCHMQPCLEGEECKVLPDRKGWS CSSGNKVKTTTRVT H
Predicted molecular weight	11 kDa
Amino acids	31 to 131

Specifications

Our **Abpromise guarantee** covers the use of **ab54415** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Lyophilized
Additional notes	Previously labelled as TAFA 2.

Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Reconstitution	For lot specific reconstitution information please contact our Scientific Support Team.

General Info

Tissue specificity	Brain-specific.
Sequence similarities	Belongs to the FAM19/TAF1 family.
Cellular localization	Cytoplasm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors