

Product datasheet

Recombinant Human FSH beta protein ab158468

1 Image

Description

Product name	Recombinant Human FSH beta protein
Expression system	Wheat germ
Accession	P01225
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	NSCELTNITAEKEEERFCISINTTWCAGYCYTRDLVYKDPA RPKIQKT CTFKELVYETVRVPGCAHHADSLYTPVATQCHCGKCDS DSTDCTVRGLG PSYCSFGEMKE
Amino acids	19 to 129
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab158468** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot
	ELISA

Form	Liquid
------	--------

Additional notes

Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
-----------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------

General Info


Relevance

Follicle stimulating hormone enables ovarian folliculogenesis to the antral follicle stage and is essential for Sertoli cell proliferation and maintenance of sperm quality in the testis. Members of the pituitary glycoprotein hormone family, of which FSH is one consist of a shared alpha chain and a beta chain encoded by a separate gene. The FSH beta gene encodes the beta subunit of follicle stimulating hormone.

Cellular localization

Secreted

Images



ab158468 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human FSH beta protein
(ab158468)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors