

Product datasheet

Recombinant Human GATAD2B protein ab153450

1 Image

Description

Product name	Recombinant Human GATAD2B protein
Expression system	Wheat germ
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Sequence	RMTE DALRLNLLKRS LDPAD ERDDVLAKRLKMEGHEAME RLKMLALLKRK DLANLEVPHELPTKQD GSGVKGYEEKLNGNLRPHGDNRT AGRPGKENIND EPVDMSAR
Amino acids	3 to 110
Tags	GST tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab153450** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	This product was previously labelled as p66 beta.

Preparation and Storage


Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

General info

Function	Has transcriptional repressor activity. Targets MBD3 to discrete loci in the nucleus.
Tissue specificity	Ubiquitous.
Sequence similarities	Contains 1 GATA-type zinc finger.
Domain	CR1 is required for interaction with MBD3 and transcription repression. CR2 is required for localization to discrete loci in the nucleus.
Cellular localization	Nucleus. Nuclear, in discrete foci.

Images



ab153450 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors