

Recombinant Human GOSR1/GS28 protein ab160530

1 Image

Description	
Product name	Recombinant Human GOSR1/GS28 protein
Expression system	Wheat germ
Accession	<u>O95249</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MAAGTSSYWEDLRKQARQLENELDKLVSFSLCTSYSH SSTRDGRRDRY SSDTTPLLNGSSQDRMFETMAIEIEQLLARLTGVNDKMAE YTNSAGVPSL NAALMHTLQRHRDILQDYTHEFHKTKANFMAIRERENLMG SVRKDIESYK SGSGVNNRRTFLKEHDHLRNSDRLEETISIAMATKENM TSQRGMLKS IHSKMNTLANRFPVNSLIQRINLRKRRDSLILGGVIGICTILL LLYAFH
Amino acids	1 to 250
Tags	GST tag N-Terminus

Specifications	
Our Abpromise guarantee covers the use of ab160530 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	Western blot ELISA
Form	Liquid
Additional notes	This product was previously labelled as GOSR1.

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info


Relevance

GOSR1/GS28 is a trafficking membrane protein which transports proteins among the endoplasmic reticulum and the Golgi and between Golgi compartments; it is considered an essential component of the Golgi SNAP receptor (SNARE) complex.

Cellular localization

Cell Membrane, Endoplasmic reticulum and Golgi Apparatus

Images



SDS-PAGE - Recombinant Human GOSR1/GS28 protein (ab160530)

ab160530 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors