

Recombinant Human GPCR C5L2 protein ab153291

1 Image

Description	
Product name	Recombinant Human GPCR C5L2 protein
Expression system	Wheat germ
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGNDSVSYEYGDYSDLSDRPVDCLDGACLAIDPLRVAPL PLYAAIFLVGV PGNAMVAWVAGKVARRRVGATWLLHLAVADLLCCLSLPI LAVPIARGGHW PYGAVGCRALPSIILLTMYASVLLLAALSADLCFLALGPAW WSTVQRACG VQVACGAAWTLALLLTVPSAIYRRLHQEHFPARLQCVVDY GGSSSTENAV TAIRFLFGFLGPLVAVASCHSALLCWAARRCRPLGTAMV GFFVCWAPYH LLGLVLTVAAPNSALLARALRAEPLVGLALAHSCLNPMFL LYFGRAQLR RSLPAACHWALRESQGQDESVDSSKKSTSHDLVSEMEV
Amino acids	1 to 337
Tags	GST tag N-Terminus

Specifications	
Our Abpromise guarantee covers the use of ab153291 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
Applications	Western blot ELISA
Form	Liquid
Additional notes	

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function

Receptor for the chemotactic and inflammatory peptide anaphylatoxin C5a, C4a and C3a and their desarginated derivatives. Couples weakly to G(i)-mediated signaling pathways.

Tissue specificity

Frontal cortex, hippocampus, hypothalamus, pons and liver.


Sequence similarities

Belongs to the G-protein coupled receptor 1 family.

Cellular localization

Cell membrane.

Images



SDS-PAGE - Recombinant Human GPCR C5L2 protein (ab153291)

ab153291 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors