

Product datasheet

Recombinant Human GTSF1 protein ab180332

1 Image

Description

Product name	Recombinant Human GTSF1 protein
Purity	> 90 % SDS-PAGE. ab180332 is purified using conventional chromatography techniques.
Expression system	Escherichia coli
Accession	<u>Q8WW33</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHH SSGLVPRGSH MGSMEETYTD SLDPEKLLQC PYDKNHQIRA CRFPYHLIKC RKNHPDVASK LATCPFNARH QVPRAEISHH ISSCDDRSCI EQDVVNQTRS LRQETLAEST WQCPCDEDW DKDLWEQTST PFVWGTTHYS DNNSPASNIV TEHKNNLASG MRVPKSLPYV LPWKNNNGNAQ
Predicted molecular weight	22 kDa including tags
Amino acids	1 to 167
Tags	His tag N-Terminus
Additional sequence information	(NP_653195).

Specifications

Our **Abpromise guarantee** covers the use of **ab180332** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Mass Spectrometry SDS-PAGE
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.32% Tris HCl, 0.88% Sodium chloride, 20% Glycerol (glycerin, glycerine), 0.02% DTT


General Info

Sequence similarities

Belongs to the UPF0224 (FAM112) family.

Contains 1 CHHC-type zinc finger.

Images



15% SDS-PAGE analysis of ab180332 (3µg).

SDS-PAGE - Recombinant Human GTSF1 protein
(ab180332)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors