

Product datasheet

Recombinant Human H-FABP protein ab52149

1 Image

Description

Product name	Recombinant Human H-FABP protein
Purity	> 95 % SDS-PAGE.
Expression system	Escherichia coli
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Amino acids	1 to 133
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab52149** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid


Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.316% Tris HCl, 50% Glycerol
------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------

General Info

Function	FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters.
Sequence similarities	Belongs to the calycin superfamily. Fatty-acid binding protein (FABP) family.
Domain	Forms a beta-barrel structure that accommodates the hydrophobic ligand in its interior.

Images



ab52149 (0.5 - 1ug/lane) on SDS-PAGE with markers.

SDS-PAGE - Recombinant Human H-FABP protein
(ab52149)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors