

Product datasheet

Recombinant Human IFNA14 protein ab124549

1 Image

Description

Product name	Recombinant Human IFNA14 protein
Purity	> 90 % SDS-PAGE. Purified by using conventional chromatography.
Expression system	Escherichia coli
Accession	<u>P01570</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSHMGSHMCNLSQTHSLNNRR TLMLMAQMRRIS PFSCCLKDRHDFEFPQEEFDGNQFQKAQAISVLHEMMQQT FNLFSTKNSSA AWDETLLEKFYIELFQQMNDLEACVIEVGVEETPLMNED SILAVKKYFQ RITLYLMEKKYSPCAWEVVRAEIMRSLSFSTNLQKRLRRK D
Predicted molecular weight	22 kDa including tags
Amino acids	24 to 189
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab124549** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.03% DTT, 0.32% Tris HCl, 50% Glycerol (glycerin, glycerine), 1.17% Sodium chloride

General Info**Function**

Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

Sequence similarities

Belongs to the alpha/beta interferon family.

Cellular localization

Secreted.

Images

15% SDS-PAGE (3 µg loaded)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors