

Recombinant Human LRPAP1 protein (FITC) ab184738

Description

Product name	Recombinant Human LRPAP1 protein (FITC)
Purity	> 90 % SDS-PAGE.
Expression system	Escherichia coli
Accession	<u>P30533</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Predicted molecular weight	39 kDa
Conjugation	FITC. Ex: 493nm, Em: 528nm

Specifications

Our **Abpromise guarantee** covers the use of **ab184738** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid
Additional notes	Extinction Coefficient = 0.93.

Preparation and Storage

Stability and Storage	Shipped on Dry Ice. Store at -80°C. Store In the Dark. pH: 6.60 Constituents: 0.58% Sodium chloride, 0.03% EDTA, 0.71% Sodium phosphate
------------------------------	---

General Info

Function	Interacts with LRP1/alpha-2-macroglobulin receptor and glycoprotein 330.
Involvement in disease	Note=In complex with the alpha-2-MR or gp330, it may have some role in the pathogenesis of membrane glomerular nephritis.

Sequence similarities	Belongs to the alpha-2-MRAP family.
Post-translational modifications	N-glycosylated.
Cellular localization	Endoplasmic reticulum. Cytoplasm. Cell surface. Intracellular and associated with cell surface receptors. Found in the endoplasmic reticulum.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors