

Product datasheet

Recombinant Human LSM4 protein ab101871

1 Image

Description

Product name	Recombinant Human LSM4 protein
Purity	> 90 % SDS-PAGE. ab101871 was purified using conventional chromatography.
Expression system	Escherichia coli
Accession	<u>Q9Y4Z0</u>
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	MGSSHHHHHHSSGLVPRGSHMLPLSLLKTAQNHPMLVE LKNGETYNGHLV SCDNWMNINLREVICTSRDGDKFWRMPECYIRGSTIKYLR PDEIIDMVK EEVVAKGRGRGGLQQKQQKGRGMGGAGRGVFGGRGR GGIPGTGRGQPEK KPGRQAGKQ
Predicted molecular weight	18 kDa including tags
Amino acids	1 to 139
Tags	His tag N-Terminus

Specifications

Our **Abpromise guarantee** covers the use of **ab101871** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE Mass Spectrometry
Mass spectrometry	MALDI-TOF
Form	Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

General Info

Function	Binds specifically to the 3'-terminal U-tract of U6 snRNA.
Sequence similarities	Belongs to the snRNP Sm proteins family.
Cellular localization	Nucleus.

Images



15% SDS PAGE analysis of 3µg ab101871.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors